

DIOSSCULPT

Diosculpt is an FDA Cleared, non-invasive 1060nm diode laser that utilizes body shaping thermotherapy technology. Diosculpt targets subcutaneous adipose tissue to achieve the local reduction of stubborn fat.



Fast Paced Growth

Body Contouring



\$7.3 Billion

Value of Global
Contouring Market in
2021

IMARC

\$11.1 Billion

Expected Value of Global
Contouring Market by
2027

IMARC





Heat or Cold?

Which is Superior?



- 1 When using cooling, only the areas that touch the applicator will be affected.

This means that a client could notice substantial difference between treated and untreated areas (lumpiness/"shelving" in the treated area).

When we use heat, it gradually dissipates at the edges of the treatment areas - resulting in a very natural look

- 2 Treatment times are 1 hour per applicator with cooling **vs 25 minutes for all four applicators with our 1060nm laser**

- 3 Results come slower with cooling (3-4 months) **vs heat with our 1060nm laser (6-12 weeks)**

HEAT with our 1060NM LASER = SUPERIOR RESULTS, FASTER!



1060nm

DioSculpt is a non-invasive 1060nm diode laser that utilizes body shaping thermotherapy technology.



Body Shaping Thermotherapy Technology



Target

DioSculpt targets subcutaneous adipose tissue to achieve the local reduction of stubborn fat.



Destroy

Unlike systems that only shrink fat cells, DioSculpt eliminates the number of fat cells to achieve the preferred body shape.

Treat Large and Small areas using 1, 2, 3 or 4 Applicators in a Session.



Arms

Back

Abdomen

Flanks

Thighs

Calves



The applicators are LARGER and MORE POWERFUL

- **More Coverage. Higher Power.**

- We deliver more energy through larger applicators. This will allow you to target a larger area in each session.
 - Using optimal power delivers superior results.
 - BE EFFICIENT.
 - **All Hands OFF Deck**
 - DioSculpt has four hands-free applicators.
- This ensures that treatments are a breeze not only for the client, but for the operator as well.

- BE FREE.

- **No Consumables**

- Unlike most aesthetic devices, there are no costly consumables with DioSculpt. We only use ultrasound gel - No disposable handpieces. No patches. No pads. There are no hidden costs for operation.

- BE PROFITABLE.

- **The Opportunities are Endless**

- DioSculpt can be used all over the body. Sculpt and define the abdomen, flanks, back, arms, thighs and calves.
 - BE THOROUGH.

- **Quick and Comfortable Sessions**

- DioSculpt sessions are quick (only 25 minutes!) and comfortable. Contact cooling protects skin and keeps clients cozy. Adipose tissue temperature is raised, without damaging surrounding tissue. There is no downtime with DioSculpt.

- BE COMFY.

- **Unbelievable Results**

- Results will continue to improve over the course of the client's DioSculpt sessions.
 - The final results can be appreciated approximately 6 weeks after the last treatment.

- BE STUNNED.

Quality Workmanship



01

Diode Laser

Built to last with a history of minimal maintenance

02

High Quality Materials

Product longevity

03

Large Applicators

4 large, hands-free applicators cover more space

04

Robust Design

Resistant to damage

Efficacy and Safety are our Priorities



01

Temperature

DioSculpt raises the treatment area to 42-47 degrees Celsius

02

Contact Cooling

Contact cooling protects skin and keeps clients comfortable

03

More Coverage

4 hands-free applicators to treat multiple areas at once

04

Safe

Safe for skin types I-VI

It's a Breeze.



Quick

Treatments are 25 minutes.
Set-up takes moments.
Clients can leave immediately after the session.

Comfortable

Controlled energy and contact cooling keep treatments cozy.



No Downtime

Activity is not restricted post-treatment. Clients can continue with their day.



Clients are Sold on the Simplicity



25 Minutes



Relax during the quick, comfortable session



Transform





DIOSSCULPT



DIOSSCULPT



DIOSSCULPT



DIOSSCULPT



DIOSSCULPT



DIOSSCULPT



DIOSSCULPT